

RNF40 antibody



Catalog Number: 114760

Product name

RNF40 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

RNF40 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, mouse liver tissue, mouse lung tissue, mouse spleen tissue. Observed molecular weight by Western-blot: 114kd

Preparation

This antibody was obtained by immunization of RNF40 recombinant protein (Accession Number: NM_001286572). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114760(RNF40 antibody) at dilution of 1:300