

RNF41 antibody



Catalog Number: 114761

Product name

RNF41 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RNF41 Rabbit Polyclonal antibody. Positive WB detected in human testis tissue. Observed molecular weight by Western-blot: 36 kDa

Preparation

This antibody was obtained by immunization of RNF41 recombinant protein (Accession Number: NM_001242826). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:114761(RNF41 antibody) at dilution of 1:500