

RNPEP antibody

Catalog Number: 114771

Product name

RNPEP antibody

IP: 1:200-1:2000

IF: 1:10-1:100


Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RNPEP Rabbit Polyclonal antibody. Positive IP detected in A549 cells. Positive WB detected in HeLa cells, A549 cells, HepG2 cells. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 72 kDa

Validations



Preparation

This antibody was obtained by immunization of RNPEP recombinant protein (Accession Number: NM_020216). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114771(RNPEP antibody) at dilution of 1:500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IF, IP


Dilutions

Recommended Dilution:

WB: 1:500-1:5000



IP Result of anti-RNPEP (IP:Catalog No:114771, 4ug; Detection:Catalog No:114771 1:500) with A549 cells lysate 2000ug.



Immunofluorescent analysis of HepG2 cells using Catalog No:114771(RNPEP Antibody) at dilution

RNPEP antibody



Catalog Number: 114771

of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)