

ROBO4-Specific antibody



Catalog Number: 114776

Product name

ROBO4-Specific antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ROBO4-Specific Rabbit Polyclonal antibody.
Positive WB detected in HeLa cells, HepG2 cells.
Observed molecular weight by Western-blot:
90-170kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_019055).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114776(ROBO4-Specific antibody) at dilution of 1:300