

RON, MST1R antibody



Catalog Number: 114785

Product name

Polyclonal

RON, MST1R antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

RON, MST1R Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue. Observed molecular weight by Western-blot: 140 kDa

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of RON, MST1R recombinant protein (Accession Number: BQ936677 NM_002447). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:114785(RON, MST1R antibody) at dilution of 1:300