

# RORC antibody



Catalog Number: 114790

---

## Product name

RORC antibody

## Immunogen

[Human RORC Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

RORC Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, Jurkat cells, mouse liver tissue. Observed molecular weight by Western-blot: 58 kDa

## Preparation

This antibody was obtained by immunization of RORC recombinant protein (Accession Number: NM\_005060). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114790(RORC antibody) at dilution of 1:200