

RORC antibody



Catalog Number: 114790

Product name

RORC antibody

Immunogen

[Human RORC Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RORC Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, Jurkat cells, mouse liver tissue. Observed molecular weight by Western-blot: 58 kDa

Preparation

This antibody was obtained by immunization of RORC recombinant protein (Accession Number: NM_005060). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114790(RORC antibody) at dilution of 1:200