

RPA1 antibody



Catalog Number: 114793

Product name

RPA1 antibody

Polyclonal

Ig Type

Rabbit IgG

Specificity

human, mouse, rat; other species not tested.

Applications

ELISA, WB

Antibody description

RPA1 Rabbit Polyclonal antibody. Positive WB detected in A431 cells, Y79 cells. Observed molecular weight by Western-blot: 70kd

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of RPA1 recombinant protein (Accession Number: NM_002945). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:114793(RPA1 antibody) at dilution of 1:200

Clonality