

RPAP1 antibody

Catalog Number: 114797

Product name

RPAP1 antibody

Immunogen

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RPAP1 Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells, HEK-293 cells, human testis tissue, mouse brain tissue, SH-SY5Y cells. Observed molecular weight by Western-blot: 153 kDa

Preparation

This antibody was obtained by immunization of RPAP1 recombinant protein (Accession Number: NM_015540). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114797(RPAP1 antibody) at dilution of 1:800



IP Result of anti-RPAP1 (IP:Catalog No:114797, 3ug; Detection:Catalog No:114797 1:800) with HeLa cells lysate 3000ug.