

# RPAP3 antibody



Catalog Number: 114799

## Product name

Polyclonal

RPAP3 antibody

## Ig Type

## Specificity

Rabbit IgG

Human; other species not tested.

## Applications

## Antibody description

ELISA, WB

RPAP3 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells, Jurkat cells. Observed molecular weight by Western-blot: 76 kDa

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Preparation

This antibody was obtained by immunization of RPAP3 recombinant protein (Accession Number: NM\_001146075). Purification method: Antigen Affinity purified.

## Validations


## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114799(RPAP3 antibody) at dilution of 1:1000