

RPN2 antibody

Catalog Number: 114820

Product name

RPN2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RPN2 Rabbit Polyclonal antibody. Positive IHC detected in human ovary tumor tissue. Positive WB detected in NIH/3T3 cells. Observed molecular weight by Western-blot: 68-69 kDa

Preparation

This antibody was obtained by immunization of RPN2 recombinant protein (Accession Number: NM_002951). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:114820(RPN2 antibody) at dilution of 1:100



(under 25x lens)

NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:114820(RPN2 Antibody) at dilution of 1:600