

# RTKN antibody



Catalog Number: 114859

## Product name

RTKN antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

RTKN Rabbit Polyclonal antibody. Positive IHC detected in human prostate cancer tissue. Positive WB detected in A431 cells, mouse brain tissue, rat brain tissue. Observed molecular weight by Western-blot: 63kd

## Preparation

This antibody was obtained by immunization of RTKN recombinant protein (Accession Number: XM\_017004635). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:114859(RTKN antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:114859(RTKN antibody) at dilution of 1:100 (under 10x lens)