

# RTN2 antibody

Catalog Number: 114862

## Product name

RTN2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

RTN2 Rabbit Polyclonal antibody. Positive IHC detected in human skeletal muscle tissue. Positive WB detected in mouse skeletal muscle tissue, human skeletal muscle tissue. Observed molecular weight by Western-blot: 20-23kd

## Preparation

This antibody was obtained by immunization of RTN2 recombinant protein (Accession Number: BC014244). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:114862(RTN2 Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human skeletal muscle slide using Catalog No:114862(RTN2 Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human skeletal muscle slide using Catalog No:114862(RTN2 Antibody) at dilution of 1:50