

RPL22 antibody

Catalog Number: 114877

Product name

RPL22 antibody

Specificity

Human; other species not tested.

Antibody description

RPL22 Rabbit Polyclonal antibody. Positive WB detected in A431 cells, HeLa cells, HepG2 cells, Jurkat cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 15-18 kDa

Preparation

This antibody was obtained by immunization of RPL22 recombinant protein (Accession Number: NM_000983). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:114877(RPL22 Antibody) at dilution of 1:600



Immunofluorescent analysis of HeLa cells using Catalog No:114877(RPL22 Antibody) at dilution of 1:50 and Rhodamine-Goat anti-Rabbit IgG