

# RPS6KA1 antibody



Catalog Number: 114912

## Product name

RPS6KA1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

RPS6KA1 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, MCF-7 cells. Observed molecular weight by Western-blot: 73 kDa

## Preparation

This antibody was obtained by immunization of RPS6KA1 recombinant protein (Accession Number: NM\_002953). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:114912(RPS6KA1 Antibody) at dilution of 1:300