

RPS6KL1 antibody



Catalog Number: 114914

Product name

RPS6KL1 antibody

Specificity

Human, hamster; other species not tested.

Antibody description

RPS6KL1 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue, human colon tissue. Positive WB detected in CHO cells. Observed molecular weight by Western-blot: 59 kDa

Preparation

This antibody was obtained by immunization of RPS6KL1 recombinant protein (Accession Number: BC004540). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



CHO cells were subjected to SDS PAGE followed by western blot with Catalog No:114914(RPS6KL1 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:114914(RPS6KL1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:114914(RPS6KL1 Antibody) at dilution of 1:50

RPS6KL1 antibody



Catalog Number: 114914

(under 40x lens)