

RRP1 antibody



Catalog Number: 114927

Product name

RRP1 antibody

Specificity

Human; other species not tested.

Antibody description

RRP1 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells, human spleen tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 53 kDa

Preparation

This antibody was obtained by immunization of RRP1 recombinant protein (Accession Number: NM_003683). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:20-1:200

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114927(RRP1 antibody) at dilution of 1:500



Immunofluorescent analysis of HeLa cells using Catalog No:114927(RRP1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)