

# RRP7A antibody



Catalog Number: 114929

## Product name

RRP7A antibody

IHC: 1:20-1:200

IF: 1:10-1:100

## Specificity

Human, Mouse; other species not tested.

## Antibody description

RRP7A Rabbit Polyclonal antibody. Positive IF detected in A431 cells. Positive IHC detected in human brain tissue, human stomach cancer tissue. Positive WB detected in mouse thymus tissue. Observed molecular weight by Western-blot: 42 kDa

## Preparation

This antibody was obtained by immunization of RRP7A recombinant protein (Accession Number: NM\_015703). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IF, IHC

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000


## Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:114929(RRP7A Antibody) at dilution of 1:1000



Immunofluorescent analysis of A431 cells using Catalog No:114929(RRP7A Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG



Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog

# RRP7A antibody



Catalog Number: 114929

---

No:114929(RRP7A Antibody) at dilution of 1:50 (under 40x lens)