

RXFP1 antibody

Catalog Number: 114943

Product name

RXFP1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RXFP1 Rabbit Polyclonal antibody. Positive IHC detected in human brain tissue, human testis tissue.

Preparation

This antibody was obtained by immunization of RXFP1 recombinant protein (Accession Number: NM_021634). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IHC


Dilutions

Recommended Dilution:

IHC: 1:50-1:500

Validations

Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:114943(RXFP1 Antibody) at dilution of 1:100 (under 10x lens).



Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:114943(RXFP1 Antibody) at dilution of 1:100 (under 40x lens).