

RXFP1 antibody



Catalog Number: 114943

Product name

RXFP1 antibody

ELISA, IHC

Specificity

Human, Mouse, Rat; other species not tested.

Dilutions


Recommended Dilution:

IHC: 1:50-1:500

Antibody description

RXFP1 Rabbit Polyclonal antibody. Positive IHC detected in human brain tissue, human testis tissue.

Validations




Preparation

This antibody was obtained by immunization of RXFP1 recombinant protein (Accession Number: NM_021634). Purification method: Antigen affinity purified.

Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:114943(RXFP1 Antibody) at dilution of 1:100 (under 10x lens).

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:114943(RXFP1 Antibody) at dilution of 1:100 (under 40x lens).

Applications