

S100A1 antibody



Catalog Number: 114953

Product name

Polyclonal

S100A1 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

S100A1 Rabbit Polyclonal antibody. Positive WB detected in human heart tissue. Observed molecular weight by Western-blot: 10 kDa

Dilutions

Recommended Dilution:

Preparation


WB: 1:200-1:2000

This antibody was obtained by immunization of S100A1 recombinant protein (Accession Number: NM_006271). Purification method: Antigen affinity purified.

Validations

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:114953(S100A1 antibody) at dilution of 1:500

Clonality