

SATB1 antibody

Catalog Number: 114977

Product name

SATB1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SATB1 Rabbit Polyclonal antibody. Positive IP detected in HEK-293 cells. Positive WB detected in human brain tissue, HEK-293 cells. Observed molecular weight by Western-blot: 98kd

Preparation

This antibody was obtained by immunization of SATB1 recombinant protein (Accession Number: NM_001322875). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:200-1:2000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114977(SATB1 antibody) at dilution of 1:4000



IP Result of anti-SATB1 (IP:Catalog No:114977, 4ug; Detection:Catalog No:114977 1:500) with HEK-293 cells lysate 2000ug.