

SATB2 antibody



Catalog Number: 114978

Product name

SATB2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SATB2 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, mouse liver tissue. Observed molecular weight by Western-blot: 80-90 kDa

Preparation

This antibody was obtained by immunization of SATB2 recombinant protein (Accession Number: NM_001172517). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114978(SATB2 antibody) at dilution of 1:1000