

# SCC4 antibody



Catalog Number: 114987

## Product name

SCC4 antibody

## Immunogen

[Human KIAA0892 Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

SCC4 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, HeLa cells. Observed molecular weight by Western-blot: 66-69 kDa

## Preparation

This antibody was obtained by immunization of SCC4 recombinant protein (Accession Number: BC063863). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:114987(KIAA0892 Antibody) at dilution of 1:600