

SCHIP1 antibody



Catalog Number: 114996

Product name

SCHIP1 antibody

Immunogen

[Human SCHIP1 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SCHIP1 Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive IHC detected in mouse brain tissue, human cerebellum tissue. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 53 kDa

Preparation

This antibody was obtained by immunization of SCHIP1 recombinant protein (Accession Number: NM_001197109). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


IHC: 1:20-1:200

IF: 1:20-1:200

Validations





HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114996(SCHIP1 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:114996(SCHIP1 Antibody) at dilution of 1:50 (under 10x lens)

SCHIP1 antibody

Catalog Number: 114996



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:114996(SCHIP1 Antibody) at dilution of 1:50 (under 40x lens)

Immunofluorescent analysis of HeLa cells using Catalog No:114996(SCHIP1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)