

SCNN1G antibody



Catalog Number: 115011

Product name

SCNN1G antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SCNN1G Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue. Positive WB detected in mouse lung tissue. Observed molecular weight by Western-blot: 70-85 kDa

Preparation

This antibody was obtained by immunization of SCNN1G recombinant protein (Accession Number: NM_001039). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:115011(SCNN1G antibody) at dilution of 1:1200



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:115011(SCNN1G antibody) at dilution of 1:100 (under 10x lens)