

# SCO1 antibody



Catalog Number: 115012

---

## Product name

SCO1 antibody

## Immunogen

## Specificity

Human; other species not tested.

## Antibody description

SCO1 Rabbit Polyclonal antibody. Positive WB detected in human heart tissue. Observed molecular weight by Western-blot: 30kd

## Preparation

This antibody was obtained by immunization of SCO1 recombinant protein (Accession Number: NM\_004589). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:115012(SCO1 antibody) at dilution of 1:400