

SCRN1 antibody



Catalog Number: 115019

Product name

SCRN1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SCRN1 Rabbit Polyclonal antibody. Positive WB detected in A549 cells, HeLa cells, mouse brain tissue. Observed molecular weight by Western-blot: 50 kDa

Preparation

This antibody was obtained by immunization of SCRN1 recombinant protein (Accession Number: NM_014766). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:115019(SCRN1 Antibody) at dilution of 1:300