

# SCRN2 antibody

Catalog Number: 115020

## Product name

SCRN2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

SCRN2 Rabbit Polyclonal antibody. Positive WB detected in BxPC-3 cells. Positive IHC detected in human gliomas tissue. Observed molecular weight by Western-blot: 46.6 kDa

## Preparation

This antibody was obtained by immunization of SCRN2 recombinant protein (Accession Number: NM\_138355). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



BxPC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:115020(SCRN2 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:115020(SCRN2 antibody) at dilution of 1:50 (under 40x lens)