

# SAFB2 antibody



Catalog Number: 115044

## Product name

SAFB2 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

SAFB2 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, BxPC-3 cells, HeLa cells, MCF7 cells, rat pancreas tissue. Positive IHC detected in human breast cancer tissue. Observed molecular weight by Western-blot: 130 kDa

## Preparation

This antibody was obtained by immunization of SAFB2 recombinant protein (Accession Number: BC025279). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115044(SAFB2 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:115044(SAFB2 antibody) at dilution of 1:50 (under 40x lens)