

SALL2 antibody



Catalog Number: 115046

Product name

SALL2 antibody

Specificity

Human; other species not tested.

Antibody description

SALL2 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HEK-293 cells. Observed molecular weight by Western-blot: 130 kDa

Preparation

This antibody was obtained by immunization of SALL2 recombinant protein (Accession Number: NM_005407). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:115046(SALL2 Antibody) at dilution of 1:1000