

# SAPK4 antibody



Catalog Number: 115061

## Product name

SAPK4 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

SAPK4 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells. Observed molecular weight by Western-blot: 42 kDa

## Preparation

This antibody was obtained by immunization of SAPK4 recombinant protein (Accession Number: NM\_002754). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:115061(MAPK13 Antibody) at dilution of 1:300