

# SDR39U1 antibody



Catalog Number: 115071

---

## Product name

SDR39U1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

SDR39U1 Rabbit Polyclonal antibody. Positive WB detected in human adrenal gland tissue. Observed molecular weight by Western-blot: 35kd

## Preparation

This antibody was obtained by immunization of SDR39U1 recombinant protein (Accession Number: NM\_020195). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



human adrenal gland tissue were subjected to SDS PAGE followed by western blot with Catalog No:115071(SDR39U1 antibody) at dilution of 1:500