

SEC14L1 antibody



Catalog Number: 115076

Product name

SEC14L1 antibody

Specificity

Human; other species not tested.

Antibody description

SEC14L1 Rabbit Polyclonal antibody. Positive WB detected in Y79 cells, HEK-293 cells. Observed molecular weight by Western-blot: 70 kDa

Preparation

This antibody was obtained by immunization of SEC14L1 recombinant protein (Accession Number: NM_001039573). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



Y79 cells were subjected to SDS PAGE followed by western blot with Catalog No:115076(SEC14L1 Antibody) at dilution of 1:300