

# SEC5/EXOC2 antibody

Catalog Number: 115082

## Product name

SEC5/EXOC2 antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

SEC5/EXOC2 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue. Positive WB detected in mouse brain tissue, human brain tissue, human ileum tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 95-100 kDa

## Preparation

This antibody was obtained by immunization of SEC5/EXOC2 recombinant protein (Accession Number: BC016918). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC, IP, WB

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000


## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115082(SEC5 antibody) at dilution of 1:500



IP Result of anti-SEC5 (IP:Catalog No:115082, 5ug; Detection:Catalog No:115082 1:500) with mouse brain tissue lysate 3000ug.



Immunohistochemical of paraffin-embedded human breast cancer using Catalog

# SEC5/EXOC2 antibody



Catalog Number: 115082

---

No:115082(SEC5 antibody) at dilution of 1:100 (under 10x lens)