

# SEC61alpha antibody



Catalog Number: 115084

## Product name

SEC61alpha antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

SEC61alpha Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in mouse brain tissue, HEK-293 cells, mouse liver tissue. Observed molecular weight by Western-blot: 52 kDa

## Preparation

This antibody was obtained by immunization of SEC61alpha recombinant protein (Accession Number: BC026179). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115084(SEC61A2 Antibody) at dilution of



Immunofluorescent analysis of HepG2 cells using Catalog No:115084(SEC61A2 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)