

# SEC63 antibody

Catalog Number: 115090

## Product name

SEC63 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

SEC63 Rabbit Polyclonal antibody. Positive WB detected in A431 cells, HEK-293 cells, PC-3 cells. Positive IP detected in A431 cells. Observed molecular weight by Western-blot: 83 kDa

## Preparation

This antibody was obtained by immunization of SEC63 recombinant protein (Accession Number: NM\_007214). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:200-1:1000

## Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:115090(SEC63 antibody) at dilution of 1:1000



IP Result of anti-SEC63 (IP:Catalog No:115090, 4ug; Detection:Catalog No:115090 1:300) with A431 cells lysate 2400ug.