

Secretin antibody



Catalog Number: 115093

Product name

Secretin antibody

Specificity

Human; other species not tested.

Antibody description

Secretin Rabbit Polyclonal antibody. Positive IHC detected in human liver cancer tissue, human stomach tissue.

Preparation

This antibody was obtained by immunization of Secretin recombinant protein (Accession Number: NM_021920). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using Catalog No:115093(SCT Antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemistry of paraffin-embedded human stomach tissue slide using Catalog No:115093(SCT Antibody) at dilution of 1:50 (under 40x lens)