

Secretin receptor antibody



Catalog Number: 115094

Product name

Secretin receptor antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Secretin receptor Rabbit Polyclonal antibody.
Positive WB detected in mouse pancreas tissue.
Observed molecular weight by Western-blot:
50-55kd

Preparation

This antibody was obtained by immunization of Secretin receptor recombinant protein (Accession Number: NM_002980). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse pancreas tissue were subjected to SDS PAGE followed by western blot with Catalog No:115094(SCTR antibody) at dilution of 1:1500