

Securin / PTTG1 antibody



Catalog Number: 115098

Product name

Securin / PTTG1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Securin / PTTG1 Rabbit Polyclonal antibody.
Positive IF detected in HepG2 cells.

Preparation

This antibody was obtained by immunization of Securin / PTTG1 recombinant protein (Accession Number: NM_001282383). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IF

Dilutions

Recommended Dilution:

IF: 1:10-1:100

Validations



Immunofluorescent analysis of HepG2 cells, using PTTG1 antibody Catalog No:115098 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).