

SELK antibody



Catalog Number: 115100

Product name

SELK antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

SELK Rabbit Polyclonal antibody. Positive WB detected in K-562 cells. Observed molecular weight by Western-blot: 10kd

Preparation

This antibody was obtained by immunization of SELK recombinant protein (Accession Number: NM_021237). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:115100(SELK antibody) at dilution of 1:400