

# SELS antibody

Catalog Number: 115102

## Product name

SELS antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

SELS Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HepG2 cells, human brain tissue. Positive IP detected in HepG2 cells. Observed molecular weight by Western-blot: 21kd

## Preparation

This antibody was obtained by immunization of SELS recombinant protein (Accession Number: NM\_018445). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115102(SELS antibody) at dilution of 1:500



IP Result of anti-SELS (IP:Catalog No:115102, 3ug; Detection:Catalog No:115102 1:500) with HepG2 cells lysate 600ug.