

# SENP3 antibody

Catalog Number: 115109

## Product name

SENP3 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

SENP3 Rabbit Polyclonal antibody. Positive WB detected in Y79 cells, HeLa cells. Positive IP detected in Y79 cells. Observed molecular weight by Western-blot: 65-70 kDa

## Preparation

This antibody was obtained by immunization of SENP3 recombinant protein (Accession Number: NM\_015670). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



Y79 cells were subjected to SDS PAGE followed by western blot with Catalog No:115109(SENP3 antibody) at dilution of 1:500



IP Result of anti-SENP3 (IP:Catalog No:115109, 4ug; Detection:Catalog No:115109 1:500) with Y79 cells lysate 2000ug.