

SEPT14 antibody



Catalog Number: 115116

Product name

SEPT14 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SEPT14 Rabbit Polyclonal antibody. Positive WB detected in rat testis tissue, mouse testis tissue. Observed molecular weight by Western-blot: 50 kDa

Preparation

This antibody was obtained by immunization of SEPT14 recombinant protein (Accession Number: BC153166). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



rat testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:115116(SEPT14 Antibody) at dilution of 1:1000