

SEPT7 antibody



Catalog Number: 115121

Product name

SEPT7 antibody

Recommended Dilution:

WB: 1:1000-1:10000

Specificity

Human, Mouse, Rat; other species not tested.


IHC: 1:20-1:200

IF: 1:20-1:200

Antibody description

SEPT7 Rabbit Polyclonal antibody. Positive IF detected in hTERT-RPE cells, HepG2 cells. Positive IHC detected in human lung cancer tissue, human gliomas tissue. Positive WB detected in HepG2 cells, A375 cells, HeLa cells, Jurkat cells, mouse testis tissue, rat testis tissue. Observed molecular weight by Western-blot: 45-50 kDa

Validations



Preparation

This antibody was obtained by immunization of SEPT7 recombinant protein (Accession Number: BC025987). Purification method: Antigen affinity purified.

HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:115121(SEPT7 antibody) at dilution of 1:2000

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal


Ig Type

Rabbit IgG

Applications

ELISA, IHC, IF, WB


Dilutions




IF result of anti-SEPT7(Catalog No:115121, 1:1000) with hTERT-RPE cell by Dr. Moshe Kim. Green: SEPT7. Blue: DNA. 1% PFA fixed

SEPT7 antibody

Catalog Number: 115121



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:115121(SEPT7 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:115121(SEPT7 antibody) at dilution of 1:100 (under 40x lens)