

SFRS2 antibody



Catalog Number: 115149

Product name

SFRS2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SFRS2 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue, A549 cells, mouse testis tissue, rat brain tissue. Observed molecular weight by Western-blot: 25-35 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_003016).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP

Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:1000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115149(SFRS2 antibody) at dilution of 1:300



IP Result of anti-SFRS2 (IP:Catalog No:115149, 4ug; Detection:Catalog No:115149 1:300) with mouse brain tissue lysate 4000ug.