

# SH3BP1 antibody



Catalog Number: 115188

## Product name

SH3BP1 antibody

## Specificity

human, mouse; other species not tested.

## Antibody description

SH3BP1 Rabbit Polyclonal antibody. Positive IF detected in HepG2. Positive WB detected in mouse lung tissue, HeLa cells, HepG2 cells, mouse brain tissue, mouse spleen tissue. Observed molecular weight by Western-blot: 76-80kd,51kd

## Preparation

This antibody was obtained by immunization of SH3BP1 recombinant protein (Accession Number: BC008282). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:10-1:100

## Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:115188(SH3BP1 antibody) at dilution of 1:500



Immunofluorescent analysis of HepG2 cells, using SH3BP1 antibody 20541-AP at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).