

SERPIND1 antibody

Catalog Number: 115207

Product name

SERPIND1 antibody

WB: 1:200-1:2000

IHC: 1:20-1:200

Immunogen

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SERPIND1 Rabbit Polyclonal antibody. Positive IHC detected in human liver cancer tissue, human hepat°Cirrhosis tissue. Positive WB detected in human plasma tissue. Observed molecular weight by Western-blot: 57-60 kDa

Preparation

This antibody was obtained by immunization of SERPIND1 recombinant protein (Accession Number: NM_000185). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:

Validations




human plasma tissue were subjected to SDS PAGE followed by western blot with Catalog No:115207(SERPIND1 Antibody) at dilution of



1:600

Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using Catalog No:115207(SERPIND1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded

SERPIND1 antibody



Catalog Number: 115207

human liver cancer tissue slide using Catalog

No:115207(SERPIND1 Antibody) at dilution of
1:50 (under 40x lens)