

SESN1 antibody



Catalog Number: 115214

Product name

SESN1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SESN1 Rabbit Polyclonal antibody. Positive WB detected in human skeletal muscle tissue, human testis tissue, mouse liver tissue. Observed molecular weight by Western-blot: 66-68 kDa

Preparation

This antibody was obtained by immunization of SESN1 recombinant protein (Accession Number: NM_014454). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:115214(SESN1 antibody) at dilution of 1:500