

# SIGLEC6 antibody

Catalog Number: 115225

## Product name

SIGLEC6 antibody

## Specificity

Human; other species not tested.

## Antibody description

SIGLEC6 Rabbit Polyclonal antibody. Positive WB detected in human serum tissue. Positive IHC detected in human placenta tissue. Observed molecular weight by Western-blot: 49kd

## Preparation

This antibody was obtained by immunization of SIGLEC6 recombinant protein (Accession Number: BC035359). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



human serum tissue were subjected to SDS PAGE followed by western blot with Catalog No:115225(SIGLEC6 antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human placenta tissue slide using Catalog No:115225(SIGLEC6 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human placenta tissue slide using Catalog No:115225(SIGLEC6 Antibody) at dilution of 1:50 (under 40x lens)